

Product Datasheet

Human Gamma-Aminobutyric Acid A Receptor Alpha 2 (gABRa2) ELISA Kit (orb1205507)

Description	Human Gamma-Aminobutyric Acid A Receptor Alpha 2 (gABRa2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P47869
Sensitivity	0.059ng/mL
Entrez	2555
Expiration Date	Please enquire.