

Product Datasheet

Human Fucosidase Alpha L2, Plasma (FUCa2) ELISA Kit (orb1205517)

Description	Human Fucosidase Alpha L2, Plasma (FUCa2) ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9BTY2
Sensitivity	0.216ng/mL
Entrez	2519
Expiration Date	Please enquire.