

Product Datasheet

Human Fas Apoptotic Inhibitory Molecule 3 (FAIM3) ELISA Kit (orb1205560)

Description	Human Fas Apoptotic Inhibitory Molecule 3 (FAIM3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O60667
Sensitivity	0.058ng/mL
Entrez	9214
Expiration Date	Please enquire.