

Product Datasheet

Human Electron Transferring Flavoprotein Dehydrogenase (ETFDH) ELISA Kit (orb1205569)

Description	Human Electron Transferring Flavoprotein Dehydrogenase (ETFDH) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q16134
Sensitivity	0.072ng/mL
Entrez	2110
Expiration Date	Please enquire.