

Product Datasheet

Human Electron Transfer Flavoprotein Alpha Polypeptide (ETFa) ELISA Kit (orb1205570)

Description	Human Electron Transfer Flavoprotein Alpha Polypeptide (ETFa) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P13804
Sensitivity	0.052ng/mL
Entrez	2108
Expiration Date	Please enquire.