

Product Datasheet

Human Endothelial Cell Adhesion Molecule (ESAM) ELISA Kit (orb1205572)

Description	Human Endothelial Cell Adhesion Molecule (ESAM) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96AP7
Sensitivity	0.117ng/mL
Entrez	90952
Expiration Date	Please enquire.