

Product Datasheet

Human Ectonucleotide Pyrophosphatase/Phosphodiesterase 1 (ENPP1) ELISA Kit (orb1205588)

Description	Human Ectonucleotide Pyrophosphatase/Phosphodiesterase 1 (ENPP1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P22413
Sensitivity	0.062ng/mL
Entrez	5167
Expiration Date	Please enquire.