

Product Datasheet

Human Elastin Microfibril Interface Located Protein 2 (EMILIN2) ELISA Kit (orb1205593)

Description	Human Elastin Microfibril Interface Located Protein 2 (EMILIN2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9BXX0
Sensitivity	0.058ng/mL
Entrez	84034
Expiration Date	Please enquire.