

Product Datasheet

Human Eukaryotic Translation Elongation Factor 1 Gamma (EEF1g) ELISA Kit (orb1205605)

Description	Human Eukaryotic Translation Elongation Factor 1 Gamma (EEF1g) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P26641
Sensitivity	0.114ng/mL
Entrez	1937
Expiration Date	Please enquire.