

Product Datasheet

Human Excitatory Amino Acid Transporter 3 (EAAT3) ELISA Kit (orb1205610)

Description	Human Excitatory Amino Acid Transporter 3 (EAAT3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P43005
Sensitivity	0.117ng/mL
Entrez	6505
Expiration Date	Please enquire.