

## Product Datasheet

### Human Excitatory Amino Acid Transporter 2 (EAAT2) ELISA Kit (orb1205611)

|                          |                                                              |
|--------------------------|--------------------------------------------------------------|
| <b>Description</b>       | Human Excitatory Amino Acid Transporter 2 (EAAT2) ELISA Kit  |
| <b>Reactivity</b>        | Human                                                        |
| <b>Range</b>             | 0.312-20ng/mL                                                |
| <b>Note</b>              | For research use only                                        |
| <b>Application notes</b> | sandwich                                                     |
| <b>Sample Types</b>      | tissue homogenates, cell lysates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>P43004</b>                                                |
| <b>Sensitivity</b>       | 0.113ng/mL                                                   |
| <b>Entrez</b>            | <b>6506</b>                                                  |
| <b>Expiration Date</b>   | Please enquire.                                              |