

Product Datasheet

Human Dolichyl Phosphate-N-Acetylglucosaminephosphotransferase 1 (DPAGT1) ELISA Kit (orb1205642)

Description	Human Dolichyl Phosphate-N-Acetylglucosaminephosphotransferase 1 (DPAGT1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H3H5
Sensitivity	0.059ng/mL
Entrez	1798
Expiration Date	Please enquire.