

## Product Datasheet

# Human Dihydroorotate Dehydrogenase (DHODH) ELISA Kit (orb1205663)

|                          |                                                                    |
|--------------------------|--------------------------------------------------------------------|
| <b>Description</b>       | Human Dihydroorotate<br>Dehydrogenase (DHODH)<br>ELISA Kit         |
| <b>Reactivity</b>        | Human                                                              |
| <b>Range</b>             | 0.156-10ng/mL                                                      |
| <b>Note</b>              | For research use only                                              |
| <b>Application notes</b> | sandwich                                                           |
| <b>Sample Types</b>      | tissue homogenates, cell<br>lysates or other biological<br>fluids. |
| <b>Uniprot ID</b>        | <b>Q02127</b>                                                      |
| <b>Sensitivity</b>       | 0.051ng/mL                                                         |
| <b>Entrez</b>            | <b>1723</b>                                                        |
| <b>Expiration Date</b>   | Please enquire.                                                    |