

Product Datasheet

Human Death Associated Protein Kinase 1 (DAPK1) ELISA Kit (orb1205684)

Description	Human Death Associated Protein Kinase 1 (DAPK1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P53355
Sensitivity	0.052ng/mL
Entrez	1612
Expiration Date	Please enquire.