

Product Datasheet

Human Dystrophin Associated Glycoprotein 1 (DAG1) ELISA Kit (orb1205687)

Description	Human Dystrophin Associated Glycoprotein 1 (DAG1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q14118
Sensitivity	0.086ng/mL
Entrez	1605
Expiration Date	Please enquire.