

Product Datasheet

Human Cysteinyl Leukotriene Receptor 1 (CYSLTR1) ELISA Kit (orb1205690)

Description	Human Cysteinyl Leukotriene Receptor 1 (CYSLTR1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9Y271
Sensitivity	0.135ng/mL
Entrez	10800
Expiration Date	Please enquire.