

Product Datasheet

Human Cytochrome P450 3A4 (CYP3A4) ELISA Kit (orb1205692)

Description	Human Cytochrome P450 3A4 (CYP3A4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P08684
Sensitivity	0.128ng/mL
Entrez	1576
Expiration Date	Please enquire.