

Product Datasheet

Human Cytochrome P450 2D6 (CYP2D6) ELISA Kit (orb1205693)

Description	Human Cytochrome P450 2D6 (CYP2D6) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P10635
Sensitivity	0.049ng/mL
Entrez	1565
Expiration Date	Please enquire.