

Product Datasheet

Human Cytochrome P450 24A1 (CYP24A1) ELISA Kit (orb1205699)

Description	Human Cytochrome P450 24A1 (CYP24A1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q07973
Sensitivity	0.053ng/mL
Entrez	1591
Expiration Date	Please enquire.