

Product Datasheet

Human Cytochrome P450 1B1 (CYP1B1) ELISA Kit (orb1205700)

Description	Human Cytochrome P450 1B1 (CYP1B1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q16678
Sensitivity	0.063ng/mL
Entrez	1545
Expiration Date	Please enquire.