

Product Datasheet

Human Cytochrome P450 1A1 (CYP1A1) ELISA Kit (orb1205702)

Description	Human Cytochrome P450 1A1 (CYP1A1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P04798
Sensitivity	0.056ng/mL
Entrez	1543
Expiration Date	Please enquire.