

Product Datasheet

Human Cysteine Rich Protein 1, Intestinal (CRIP1) ELISA Kit (orb1205741)

Description	Human Cysteine Rich Protein 1, Intestinal (CRIP1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P50238
Sensitivity	0.063ng/mL
Entrez	1396
Expiration Date	Please enquire.