

Product Datasheet

Human Cellular Retinoic Acid Binding Protein 2 (CRABP2) ELISA Kit (orb1205750)

Description	Human Cellular Retinoic Acid Binding Protein 2 (CRABP2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P29373
Sensitivity	0.114ng/mL
Entrez	1382
Expiration Date	Please enquire.