

## Product Datasheet

### Human Carboxypeptidase A3, Mast Cell (CPA3) ELISA Kit (orb1205764)

|                          |                                                          |
|--------------------------|----------------------------------------------------------|
| <b>Description</b>       | Human Carboxypeptidase A3,<br>Mast Cell (CPA3) ELISA Kit |
| <b>Reactivity</b>        | Human                                                    |
| <b>Range</b>             | 78.1-5000pg/mL                                           |
| <b>Note</b>              | For research use only                                    |
| <b>Application notes</b> | sandwich                                                 |
| <b>Sample Types</b>      | serum, plasma or other<br>biological fluids.             |
| <b>Uniprot ID</b>        | <b>P15088</b>                                            |
| <b>Sensitivity</b>       | 34pg/mL                                                  |
| <b>Entrez</b>            | <b>1359</b>                                              |
| <b>Expiration Date</b>   | Please enquire.                                          |