

Product Datasheet

Human Collagen Type VIII Alpha 1 (COL8a1) ELISA Kit (orb1205777)

Description	Human Collagen Type VIII Alpha 1 (COL8a1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P27658
Sensitivity	0.052ng/mL
Entrez	1295
Expiration Date	Please enquire.