

Product Datasheet

Human Collagen Type VI Alpha 3 (COL6a3) ELISA Kit (orb1205779)

Description	Human Collagen Type VI Alpha 3 (COL6a3) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P12111
Sensitivity	12.7pg/mL
Entrez	1293
Expiration Date	Please enquire.