

Product Datasheet

Human Collagen Type IV Alpha 3 (COL4a3) ELISA Kit (orb1205782)

Description	Human Collagen Type IV Alpha 3 (COL4a3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q01955
Sensitivity	0.113ng/mL
Entrez	1285
Expiration Date	Please enquire.