

Product Datasheet

Human Collagen Type III Alpha 1 (COL3a1) ELISA Kit (orb1205783)

Description	Human Collagen Type III Alpha 1 (COL3a1) ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P02461
Sensitivity	0.226ng/mL
Entrez	1281
Expiration Date	Please enquire.