

Product Datasheet

Human Chymase 1, Mast Cell (CMA1) ELISA Kit (orb1205798)

Description	Human Chymase 1, Mast Cell (CMA1) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P23946
Sensitivity	6.5pg/mL
Entrez	1215
Expiration Date	Please enquire.