

Product Datasheet

Human Chloride Intracellular Channel Protein 4 (CLIC4) ELISA Kit (orb1205800)

Description	Human Chloride Intracellular Channel Protein 4 (CLIC4) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9Y696
Sensitivity	30pg/mL
Entrez	25932
Expiration Date	Please enquire.