

Product Datasheet

Human Chloride Intracellular Channel Protein 1 (CLIC1) ELISA Kit (orb1205801)

Description	Human Chloride Intracellular Channel Protein 1 (CLIC1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	000299
Sensitivity	0.123ng/mL
Entrez	1192
Expiration Date	Please enquire.