

Product Datasheet

Human C-Type Lectin Domain Family 4, Member K (CLEC4K) ELISA Kit (orb1205802)

Description	Human C-Type Lectin Domain Family 4, Member K (CLEC4K) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9UJ71
Sensitivity	0.039ng/mL
Entrez	50489
Expiration Date	Please enquire.