

Product Datasheet

Human Charcot Leyden Crystal Protein (CLC) ELISA Kit (orb1205808)

Description	Human Charcot Leyden Crystal Protein (CLC) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q05315
Sensitivity	0.128ng/mL
Entrez	1178
Expiration Date	Please enquire.