

Product Datasheet

Human Cardiotrophin Like Cytokine Factor 1 (CLCF1) ELISA Kit (orb1205809)

Description	Human Cardiotrophin Like Cytokine Factor 1 (CLCF1) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9UBD9
Sensitivity	13.9pg/mL
Entrez	23529
Expiration Date	Please enquire.