

Product Datasheet

Human Cholinergic Receptor, Nicotinic, Alpha 7 (CHRNA7) ELISA Kit (orb1205817)

Description	Human Cholinergic Receptor, Nicotinic, Alpha 7 (CHRNA7) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P36544
Sensitivity	0.112ng/mL
Entrez	1139
Expiration Date	Please enquire.