

Product Datasheet

Human Cholinergic Receptor, Nicotinic, Alpha 3 (CHRNA3) ELISA Kit (orb1205819)

Description	Human Cholinergic Receptor, Nicotinic, Alpha 3 (CHRNA3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P32297
Sensitivity	0.107ng/mL
Entrez	1136
Expiration Date	Please enquire.