

Product Datasheet

Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 6 (CEACAM6) ELISA Kit (orb1205833)

Description	Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 6 (CEACAM6) ELISA Kit
Reactivity	Human
Range	23.4-1500pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P40199
Sensitivity	9.68pg/mL
Entrez	4680
Expiration Date	Please enquire.