

Product Datasheet

Human Cerebral Dopamine Neurotrophic Factor (CDNF) ELISA Kit (orb1205837)

Description	Human Cerebral Dopamine Neurotrophic Factor (CDNF) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q49AH0
Sensitivity	0.057ng/mL
Entrez	441549
Expiration Date	Please enquire.