

Product Datasheet

Human Cyclin Dependent Kinase Inhibitor 3 (CDKN3) ELISA Kit (orb1205838)

Description	Human Cyclin Dependent Kinase Inhibitor 3 (CDKN3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q16667
Sensitivity	0.053ng/mL
Entrez	1033
Expiration Date	Please enquire.