

Product Datasheet

Human Cadherin 17 (CDH17) ELISA Kit (orb1205842)

Description	Human Cadherin 17 (CDH17) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q12864
Sensitivity	0.054ng/mL
Entrez	1015
Expiration Date	Please enquire.