

Product Datasheet

Human Chemokine C-C-Motif Receptor 6 (CCR6) ELISA Kit (orb1205846)

Description	Human Chemokine C-C-Motif Receptor 6 (CCR6) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P51684
Sensitivity	0.053ng/mL
Entrez	1235
Expiration Date	Please enquire.