

Product Datasheet

Human Collagen and calcium-binding EGF domain-containing protein 1 (CCBE1) ELISA Kit (orb1205852)

Description	Human Collagen and calcium-binding EGF domain-containing protein 1 (CCBE1) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6UXH8
Sensitivity	11.73pg/mL
Entrez	147372
Expiration Date	Please enquire.