

Product Datasheet

Human Carbonic Anhydrase XII (CA12) ELISA Kit (orb1205882)

Description	Human Carbonic Anhydrase XII (CA12) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O43570
Sensitivity	0.124ng/mL
Entrez	771
Expiration Date	Please enquire.