

## Product Datasheet

# Human Complement Component 1, Q Subcomponent C (C1qC) ELISA Kit (orb1205885)

|                          |                                                                 |
|--------------------------|-----------------------------------------------------------------|
| <b>Description</b>       | Human Complement Component 1, Q Subcomponent C (C1qC) ELISA Kit |
| <b>Reactivity</b>        | Human                                                           |
| <b>Range</b>             | 31.2-2000pg/mL                                                  |
| <b>Note</b>              | For research use only                                           |
| <b>Application notes</b> | sandwich                                                        |
| <b>Sample Types</b>      | serum, plasma or other biological fluids.                       |
| <b>Uniprot ID</b>        | <b>P02747</b>                                                   |
| <b>Sensitivity</b>       | 14.8pg/mL                                                       |
| <b>Entrez</b>            | <b>714</b>                                                      |
| <b>Expiration Date</b>   | Please enquire.                                                 |