

Product Datasheet

Human Breast Cancer Anti-Estrogen Resistance 1 (BCAR1) ELISA Kit (orb1205905)

Description	Human Breast Cancer Anti-Estrogen Resistance 1 (BCAR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P56945
Sensitivity	0.053ng/mL
Entrez	9564
Expiration Date	Please enquire.