

Product Datasheet

Human ATPase, H⁺ Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) ELISA Kit (orb1205918)

Description	Human ATPase, H ⁺ Transporting, Lysosomal Accessory Protein 2 (ATP6AP2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	075787
Sensitivity	0.112ng/mL
Entrez	10159
Expiration Date	Please enquire.