

Product Datasheet

Human Calcium-Transporting ATPase Type 2C Member 1 (ATP2C1) ELISA Kit (orb1205920)

Description	Human Calcium-Transporting ATPase Type 2C Member 1 (ATP2C1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P98194
Sensitivity	0.077ng/mL
Entrez	27032
Expiration Date	Please enquire.