

Product Datasheet

Human ATPase, Na⁺/K⁺ Transporting Beta 1 Polypeptide (ATP1b1) ELISA Kit (orb1205923)

Description	Human ATPase, Na ⁺ /K ⁺ Transporting Beta 1 Polypeptide (ATP1b1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P05026
Sensitivity	0.066ng/mL
Entrez	481
Expiration Date	Please enquire.