

Product Datasheet

Human ATPase, Na⁺/K⁺ Transporting Alpha 2 Polypeptide (ATP1a2) ELISA Kit (orb1205924)

Description	Human ATPase, Na ⁺ /K ⁺ Transporting Alpha 2 Polypeptide (ATP1a2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P50993
Sensitivity	0.105ng/mL
Entrez	477
Expiration Date	Please enquire.