

Product Datasheet

Human Astacin Like Metallo Endopeptidase (ASTL) ELISA Kit (orb1205928)

Description	Human Astacin Like Metallo Endopeptidase (ASTL) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q6HA08
Sensitivity	0.060ng/mL
Entrez	431705
Expiration Date	Please enquire.